

# TCM5 Connected devices

**Updated:  
February 2022**

Manufacturer	Device	Connectivity available TCM5	Manufacturer SW version	From TCM5 SW	Comments	
Patient Monitoring Systems	Dräger	Infinity Delta	✓	VF10.1-W	1.3.0	
		Infinity Delta XL	✓	VF10.1-W	1.3.0	
		Infinity Kappa	✓	VF10.1-W	1.3.0	
		Infinity Omega	✓	VF10.1-W	1.3.0	
	Mindray	Beneview T5 / T6 / T8 / T9	✓		1.2.0	
	Philips	IntelliView MP40-90	✓		1.1.8	Using Philips EC10 Interface Module or IntelliBridge EC10 integral interface board with Open Interface Driver (ED101) and EC5 ID Module #102
		IntelliView MX400-800	✓		1.1.8	Using Philips EC10 Interface Module or IntelliBridge EC10 integral interface board with Open Interface Driver (ED101) and EC5 ID Module #102
	GE Healthcare	Carescape B450	✓	3.1.353 6ICU and 2 6ORP	1.3.0	
		Carescape B650	✓	3.1.353 6ICU and 2 6ORP	1.3.0	
		Carescape B850	✓	3.1.353 6ICU and 2 6ORP	1.3.0	
Dash 5000		✓	7.3 MULTI 26MAR15	1.3.0		
Nihon Kohden	BSM series	✓	V.01-03	1.2.0		
	CSM series	✓	Ver.00-08 12	1.4		
Polysomnographs	Cidelec	CID analog	✓		1.2.0	
	NOX Medical	A1	✓	Noxturnal software version 5.1.0.19071	1.4.0	
		C1	✓	Noxturnal software version 5.1.0.19071	1.4.0	
	Compumedics	E-Series	✓		1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
		Grael	✓		1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
		Siesta	✓		1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
		Somté	✓		1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later

	Manufacturer	Device	Connectivity available	Manufacturer SW version	From TCM5 SW version	Comments
Polysomnographs	Natus	Embletta Gold	✓		1.2.0	
		Embla S7000 / N7000	✓		1.2.0	
		Embla S4000	✓		1.2.0	
		Embla S4500	✓		1.2.0	
		Mdrive	✓		1.2.0	
	Philips Respironics	Alice analog	✓		1.2.0	
Ventilators	Breas	Nippy4+	✓		1.3.0	
		Vivo45	✓		1.3.0	
		Vivo45LS	✓		1.3.0	
EMR	Capsule Technologies	Medical Device Information Platform	✓		1.3.0	
CPET	Cosmed	Quark RMR	✓		1.4.0	

✓ Connection validated